

**Erudit Basic Mono 40x-400x Mikroskop
Erudit Basic Mono 40x-400x Microscope**

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ABOUT THESE INSTRUCTIONS

These operating instructions are to be regarded as part of the device! Read the safety instructions and the instruction manual carefully before using this device. Keep this instruction manual in a safe place for future reference. When the device is sold or given to someone else, the instruction manual must be provided to the new owner/user of the product.

INTENDED USE

This product is for private use only. It was developed for magnifying objects of nature.

This device is not intended to be used by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and/or knowledge, unless they have been given supervision or instruction on how to use the device by a person responsible for their safety.

Children should be supervised to ensure that they do not play with the device.

GENERAL SAFETY INFORMATION

CHOKING HAZARD!

Keep packaging materials (plastic bags, rubber bands, etc.) away from children! CHOKING HAZARD!

RISK OF INJURY!

Includes sharp-edged corners and points!

This device often requires the use of sharp-edged and pointed tools. For this reason, keep this device together with all accessories and tools in a place inaccessible to children. RISK OF INJURY!

RISK OF ELECTRIC SHOCK!

This unit contains electronic parts that are powered by a power source (batteries). Never leave children unattended when handling the device! Only use the device as described in the instructions, otherwise there is a RISK of an ELECTRIC SHOCK!

Check the device, cables and connections for damage before use.

Do not use the device if it is damaged or contains damaged live parts! Damaged parts must be replaced immediately by an authorized service centre.

FIRE/EXPLOSION HAZARD!

Do not expose the unit to high temperatures. Only use the recommended batteries. Do not short-circuit the device or batteries or throw them into a fire! Excessive heat and improper handling can cause a short circuit, fire or even an explosion!

! RISK OF PROPERTY DAMAGE!

Do not subject the device to shocks!

The manufacturer accepts no liability for damage resulting from incorrectly inserted batteries!

Do not disassemble the device! In case of a defect, contact your dealer. The dealer will contact the service centre and send the device for repair if necessary.

1. PACKAGE CONTENTS (FIG. 1)

(Some accessories are already pre-assembled)

- (A) Microscope with
 - (B) Microscope head
 - (C) 10x WF eyepiece
 - (D) Objective lenses: 4x, 10x, 40x
 - (E) Condenser
- (F) Box with 5 permanent specimens, empty slides and coverslips
- (G) Smartphone holder
- (H) Carry case

2. PARTS OVERVIEW FOR ALL MODELS (FIG. 2)

- 1 Eyepiece
- 2 Eyepiece tubes
- 3 Microscope head
- 4 Nosepiece
- 5 Objective lens (4x, 10x or 40x)
- 6 Specimen clamp
- 7 Condenser
- 8 Adjusting lever for iris diaphragm
- 9 Filter holder
- 10 Condenser height adjustment
- 11 Condenser adjustment screw
- 12 Light
- 13 Coarse focus
- 14 Fine focus
- 15 X-axis adjustment for XY stage
- 16 Y-axis adjustment for XY stage
- 17 Rotary switch (on/off switch with brightness control)
- 18 Battery compartment
- 19 Screw plug for battery compartment
- 20 Smartphone holder

only with Erudit Basic Bino (Fig. 3)

- 21 Eyepiece tubes with foldable central axis
- 22 Vision correction (dioptre adjustment)

3. LOCATION

Before setting up your microscope, choose a suitable location.

First of all, you should make sure that your microscope is placed on a stable, vibration-free surface.

4. INSERT THE BATTERIES* (FIG. 4)

The battery compartment (18) is located on the underside of the microscope. Loosen the screw plug (19) on the battery compartment cover. Remove the cover. Insert three AA-LR6 batteries*. Make sure the polarity is correct as indicated on the inside of the battery compartment. Close the battery compartment again and tighten the screw (19) to a hand-tight position only to avoid damaging the cover.

Alternatively, you can also use rechargeable batteries of the same type. However, rechargeable batteries have a lower voltage, and therefore the light may be less bright.

5. ELECTRIC LIGHT (FIG. 5)

The microscope is equipped with a battery-powered LED light. This allows you to operate the microscope without a mains connection (e.g. outdoors).

The on/off switch and brightness control are combined in one rotary switch (17). To switch on the microscope, turn the rotary switch (17) until you hear a click. To adjust the brightness, continue turning the rotary switch. Adjust the brightness of the light to obtain the best possible illumination.

! NOTE!

Do not operate the microscope at maximum brightness for a long period of time. This will reduce the lifespan of the LED.

The battery will power the light for up to 30 hours (at the lowest brightness setting).

6. OBSERVATION (FIG. 6)

Start each observation with the lowest magnification (10x eyepiece and 4x objective). This makes it easier to centre and adjust the specimen.

The higher the magnification, the more light is needed to obtain a good quality image. First, set the nosepiece (4) to the lowest magnification and insert the 10x eyepiece. Now place a specimen directly under the objective on the microscope stage and fix it in place with the specimen clamp (6). Move the specimen using the X-axis (15) and Y-axis adjustment (16) for the XY stage. Now look through the eyepiece (1) and carefully turn the coarse focus adjustment (13) until the image becomes visible. Use the fine focus (14) to obtain a sharper image.

! NOTE!

Never turn the fine focus adjustment beyond the stop position.

7. ADJUSTING THE INTERPUPILLARY DISTANCE

(Erudit Basic Bino only)

Place the specimen on the specimen stage and bring the image into focus.

Set the correct interpupillary distance for binocular observation by folding the eyepiece tubes (21) on the microscope head (3) so that the right and left fields of view merge into one when viewing.

8. CONFIGURING THE VISION CORRECTION

(Erudit Basic Bino only)

Place the specimen on the specimen stage. Turn the 40x objective (5) to the working position. First, observe only through the right eyepiece with your right eye, and then use the coarse and fine focus knobs to focus the image. Next, observe only through the left eyepiece with your left eye, and then use the dioptre adjustment ring (22) to focus the image.

9. CONDENSER

The condenser (7) concentrates the light from the LED unit onto the specimen. The condenser can be adjusted using the condenser height adjustment (10).

Precise adjustments can be made using the two adjustment screws (11). A fine slotted screwdriver or Allen key (depending on the model) is required for this. To adjust the condenser, close the iris diaphragm

(8) completely and focus on a specimen with the 4x objective (5). Turn the condenser (7) slightly downwards via the condenser height adjustment (10) until you can clearly see the edge of the aperture. If the bright area is not in the centre, adjust the condenser (7) accordingly with the two adjustment screws (11). Once this process is complete, bring the condenser (7) back into the optimal observation position and open the aperture (8) again.

10. IRIS DIAPHRAGM

The iris diaphragm (8) allows you to increase the depth of field. An open aperture lets more light through, with a shallow depth of field. A partially closed aperture increases the depth of field, but this requires increasing the brightness of the light.

The iris diaphragm (8) should not be closed too far, especially at higher magnifications, as this will reduce the resolution.

11. FILTER HOLDER

Filters (available separately) with a diameter of 32 mm and a thickness of approx. 1.5 mm can be inserted into the filter holder (9).

12. SMARTPHONE HOLDER (FIG. 7)

The smartphone holder (20) can be placed on the eyepiece and fixed in place with the knurled screw.

For the binocular model (5102200), the left or right eyepiece can be used for attaching the holder.

! NOTE!

The suction cups must be clean and free from dust and dirt. It can be helpful to slightly moisten the suction cups. Make sure that the smartphone holder is fixed securely in place and that the smartphone is firmly seated on the holding plate.

Smartphones with a rough surface hold less well than those with a smooth surface. Make sure that the smartphone does not slip off the holding plate.

Now start the camera app. The camera lens of the smartphone must be exactly above the eyepiece so that the image is exactly centred on the display.

It may be necessary to use the zoom function to fill the screen. There may be slight shading at the edges.

Remove the smartphone from the holder after use.

13. MAINTENANCE AND STORAGE

Only use the recommended batteries.

Always replace weak or used batteries with a complete new set of batteries with full capacity. Do not use batteries of different brands, types or different capacities. Remove the batteries from the device if you do not plan to use it for an extended period!

Do not disassemble the microscope!

The microscope is a precision optical instrument, and therefore must not be disassembled. This can seriously damage its functionality. Do not disassemble the objective lenses!

Avoid dirtying the lenses!

Prevent dust or moisture from coming into contact

with your microscope. Avoid fingerprints on all optical surfaces! If dust or dirt does get on your microscope or accessories, remove it as described below.

Cleaning the non-optical components

Before cleaning the non-optical components (e.g. housing, specimen stage), remove the batteries from the device. Dust or dirt on your microscope or accessories should first be removed with a soft brush. Then clean the affected area with a soft, lint-free cloth, which can be dampened with a mild detergent solution. Allow any residual moisture to dry completely before using the microscope. DO NOT use organic solvents (such as alcohol, ether, acetone, xylene or other thinners) to clean painted or plastic parts!

Cleaning the optical components

The eyepiece and objective lenses of the microscope come with a premium coating to increase the quality of the image. These parts should not be wiped dry, as dry dirt or dust can scratch the sensitive coating. Ideally, the optical elements should be removed from the microscope before cleaning. First, blow any loose dust off the lens surfaces. Then wipe the lens surface clean using a good quality lens cloth or soft cloth slightly moistened with pure alcohol (available at pharmacies).

When not in use

After use, return the microscope to the case and place it in a dry and mould-free location. If the microscope and accessories are not going to be used for a long period of time, they should be stored in their respective cases. We recommend storing all objectives and eyepieces in a closed container with desiccant. Remove the batteries from the device if you do not plan to use it for an extended period!

Remember:

A well-maintained microscope will retain its optical quality (and thus its value) for many years.

14. USEFUL FACTS

Aperture:

The numerical aperture (unit of measurement for the resolving power of the objectives) is engraved next to the magnification number on the objectives.

Calculating the microscope magnification level:

Formula:

Eyeiece magnification x objective magnification = total magnification

Example:

10x (eyepiece) x 10x (objective) = 100x total magnification

15. DISPOSAL

 Dispose of the packaging materials according to type. Observe the current legal regulations when disposing of the device! Information on proper disposal can be obtained from your municipal waste disposal service or environmental agency.

 Do not throw electrical appliances in household waste!

According to the European Directive 2002/96/EC on waste electrical and electronic equipment and its transposition into national law, used electrical devices must be collected separately and recycled in an environmentally friendly manner.

16. DECLARATION OF CONFORMITY

16.1 CE DECLARATION OF CONFORMITY

Product type: Transmitted light microscope

Product name: Erudit Basic Mono/Bino

Item no.: 5102100 / 5102200

Bresser GmbH has prepared a "Declaration of Conformity" in accordance with the applicable guidelines and corresponding standards. This can be viewed at:

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Suite 3G, Eden House · Enterprise Way
Edenbridge, Kent TN8 6HF · Great Britain

17. WARRANTY & SERVICE

The regular warranty period is 2 years, which begins on the day of purchase. To benefit from an extended voluntary warranty period as indicated on the gift box, registration on our website is required. The complete warranty conditions as well as information on the warranty extension and services can be found at www.bresser.de/warranty_terms

18. TECHNICAL DATA

	Erudit Basic Mono (Art. No. 5102100)	Erudit Basic Bino (Art. No. 5102200)
Microscope head	Monocular	Binocular
XY stage	with vernier adjustment	
Eyepiece(s)	1 pc. DIN WF 10x	1 pair DIN WF 10x
Objective lens	4x / 10x / 40x	
Magnification level	40x / 100x / 400x	
Power supply	4.5 V (3 batteries, type AA/LR6) (not included)	
Illumination	LED* (battery-operated)	

* The LED light installed in this unit has a very long lifespan and does not need to be replaced. A replacement is also not possible for technical reasons!

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BRESSER Benelux B.V.
Bresser UK Ltd.