




Sensor-/Bildfeldvergleich

Maßstab 4:1, verwendetes Mikroskop: BRESSER Science Infinity, 10x Objektiv
Objekt: Pflanzendünnschnitt


WF 10× / 22 mm Okular




Art.-Nr. 11798




Art.-Nr. 11979




Art.-Nr. 12298




Art.-Nr. 14268




Art.-Nr. 11967




Art.-Nr. 12752




Art.-Nr. 14055




Art.-Nr. 14099




Art.-Nr. 12207



Art.-Nr. 13414




Art.-Nr. 16304



Art.-Nr. 12520



Art.-Nr. 10983



Art.-Nr. 11297



Art.-Nr. 11798

1/2.7"
1920 × 1080 Pixel

1:1

Art.-Nr. 11979

1/2.8"
2048 × 1536 Pixel

1:1

Art.-Nr. 12298

1/2.5"
2529 × 1944 Pixel

1:1

Art.-Nr. 14268

1"
4096 × 2160 Pixel

1:1

Art.-Nr. 11967

2/3"
2448 × 2048 Pixel

1:1

Art.-Nr. 12752

1/1.7"
4000 × 3000 Pixel

1:1

Art.-Nr. 14055

1"
5440 × 3648 Pixel

1:1

Art.-Nr. 14099

1/2.7"
1920 × 1080 Pixel

1:1

Art.-Nr. 12207

1/1.8"
1920 × 1080 Pixel - Live
2592 × 1944 Pixel- Foto

1:1

Art.-Nr. 13414

1/1.9"
1920 × 1080 Pixel

1:1

NEU

Art.-Nr. 16304

1/1.8"
3840 × 2160 Pixel

1:1

NEU

Art.-Nr. 12520

1/2.8"
1920 × 1080 Pixel

1:1

NEU

Art.-Nr. 10983

1/2.8"
2592 × 1944 Pixel

1:1

NEU

Art.-Nr. 11297

1/0.25"
2592 × 1944 Pixel

1:1

Übersicht: MikroCams



Technische Daten

Art.-Nr.	Artikelbezeichnung	Auflösung Pixel	Pixel- größe µm	Zoll "	Länge mm	Breite mm	Fläche mm²	Diagonale mm	Seiten- verhältnis	Modus
11798	MikroOkular Full HD Okularkamera	1920 x 1080	3,0	1/2.7	5,76	3,24	18,66	6,61	16:9	
11979	MikroCam SP 3.1 Mikroskopkamera	2048 x 1536	2,5	1/2.8	5,28	3,84	20,28	6,53	4:3	
12298	MikroCam SP 5.0 Mikroskopkamera	2592 x 1944	2,2	1/2.5	5,70	4,30	24,51	7,14	4:3	
14268	MikroCamII 9MP 4K 1" Mikroskopkamera	4096 x 2160	3,45	1.0	14,13	7,45	105,27	15,97	~ 17:9	
11967	MikroCamII 5MP HIS Mikroskopkamera	2448 x 2048	3,45	2/3.0	8,45	7,07	59,74	11,02	5:4	
12752	MikroCamII Mikroskopkamera 12MP	4000 x 3000	1,85	1/1.7	7,40	5,55	41,07	9,25	4:3	
14055	MikroCamII 20 MP 1" Mikroskopkamera	5440 x 3648	2,4	1.0	13,06	8,76	114,41	15,73	3:2	
14099	MikroCam MINI Full HD HDMI Mikroskopkamera	1920 x 1080	2,9	1/2.8	5,57	3,13	17,43	6,39	16:9	
12207	MikroCam PRO HDMI Mikroskopkamera 5MP	1920 x 1080 2592 x 1944	2,4	1/1.8	6,22	4,67	29,05	7,78	16:9 4:3	Live Photo
13414	MikroCam PRO HDMI Autofocus Mikroskopkamera	1920 x 1080	3,75	1/1.9	7,20	4,05	29,16	8,26	16:9	
16304	MikroCam PRO HDMI 4K Mikroskopkamera	3840 x 2160	2,0	1/1.8	7,68	4,32	33,18	8,81	16:9	
12520	MikroCam SF Full-HD-Mikroskopkamera	1920 x 1080	2,9	1/2.8	5,57	3,13	17,43	6,39	16:9	
10983	MikroCam SF 5.0 MP Mikroskopkamera	2592 x 1944	2,0	1/2.8	5,18	3,89	20,15	6,48	4:3	
11297	MikroCam SF UHD Mikroskopkamera	3840 x 2160	1,62	1/2.5	6,22	3,50	21,77	7,14	16:9	

NEU