



Sensor-/Bildfeldvergleich

Maßstab 4:1, verwendetes Mikroskop: BRESSER Science Infinity, 10x Objektiv
Objekt: Pflanzendünnschnitt


WF 10× / 22 mm Okular




Art.-Nr. 5913650




Art.-Nr. 5914321




Art.-Nr. 5914520




Art.-Nr. 5914090




Art.-Nr. 5914510




Art.-Nr. 5914912




Art.-Nr. 5914200




Art.-Nr. 5914170




Art.-Nr. 5914185




Art.-Nr. 5914190



Art.-Nr. 5914195



Art.-Nr. 5914122



Art.-Nr. 5914522



Art.-Nr. 5914822



Art.-Nr. 5913650

1/2.7"
1920 × 1080 Pixel

1:1

Art.-Nr. 5914321

1/2.8"
2048 × 1536 Pixel

1:1

Art.-Nr. 5914520

1/2.5"
2529 × 1944 Pixel

1:1

Art.-Nr. 5914090

1"
4096 × 2160 Pixel

1:1

Art.-Nr. 5914510

2/3"
2448 × 2048 Pixel

1:1

Art.-Nr. 5914912

1/1.7"
4000 × 3000 Pixel

1:1

Art.-Nr. 5914200

1"
5440 × 3648 Pixel

1:1

Art.-Nr. 5914170

1/2.7"
1920 × 1080 Pixel

1:1

Art.-Nr. 5914185

1/1.8"
1920 × 1080 Pixel - Live
2592 × 1944 Pixel- Foto

1:1

Art.-Nr. 5914190

1/1.9"
1920 × 1080 Pixel

1:1

NEU

Art.-Nr. 5914195

1/1."
3840 × 2160 Pixel

1:1

NEU

Art.-Nr. 5914122

1/2.8"
1920 × 1080 Pixel

1:1

NEU

Art.-Nr. 5914522

1/2.8"
2592 × 1944 Pixel

1:1

NEU

Art.-Nr. 5914822

1/0.25"
2592 × 1944 Pixel

1:1

NEU

Art.-Nr. 5914250

1/2.3"
4928 × 4928 Pixel

1:1

Übersicht: MikroCams



Technische Daten

Art.-Nr.	Artikelbezeichnung	Auflösung Pixel	Pixel- größe µm	Zoll "	Länge mm	Breite mm	Fläche mm²	Diagonale mm	Seiten- verhältnis	Mo- dus
5913650	MikroOkular Full HD Okularkamera	1920 x 1080	3,0	1/2.7	5,76	3,24	18,66	6,61	16:9	
5914321	MikroCam SP 3.1 Mikroskopkamera	2048 x 1536	2,5	1/2.8	5,28	3,84	20,28	6,53	4:3	
5914520	MikroCam SP 5.0 Mikroskopkamera	2592 x 1944	2,2	1/2.5	5,70	4,30	24,51	7,14	4:3	
5914090	MikroCamII 9MP 4K 1" Mikroskopkamera	4096 x 2160	3,45	1.0	14,13	7,45	105,27	15,97	~ 17:9	
5914510	MikroCamII 5MP HIS Mikroskopkamera	2448 x 2048	3,45	2/3.0	8,45	7,07	59,74	11,02	5:4	
5914912	MikroCamII Mikroskopkamera 12MP USB 3.0	4000 x 3000	1,85	1/1.7	7,40	5,55	41,07	9,25	4:3	
5914200	MikroCamII 20 MP 1" Mikroskopkamera	5440 x 3648	2,4	1.0	13,06	8,76	114,41	15,73	3:2	
5914170	MikroCam MINI Full HD HDMI Mikroskopkamera	1920 x 1080	2,9	1/2.8	5,57	3,13	17,43	6,39	16:9	
5914185	MikroCam PRO HDMI Mikroskopkamera 5MP	1920 x 1080 2592 x 1944	2,4	1/1.8	6,22	4,67	29,05	7,78	16:9 4:3	Live Photo
5914190	MikroCam PRO HDMI Autofocus Mikroskopkamera	1920 x 1080	3,75	1/1.9	7,20	4,05	29,16	8,26	16:9	
5914195	MikroCam PRO HDMI 4K	3840 x 2160	2,0	1/1.8	7,68	4,32	33,18	8,81	16:9	
5914122	MikroCam SF Full HD	1920 x 1080	2,9	1/2.8	5,57	3,13	17,43	6,39	16:9	
5914522	MikroCam SF 5.0 MP	2592 x 1944	2,0	1/2.8	5,18	3,89	20,15	6,48	4:3	
5914822	MikroCam SF UHD	3840 x 2160	1,62	1/2.5	6,22	3,50	21,77	7,14	16:9	

NEU