



Sensor-/Bildfeldvergleich

Scale 4:1, used microscope: BRESSER Science Infinity, 10x objective lens
Object: thinning of plants


WF 10× / 22 mm Okular





Art No. 5913650




Art No. 5914321




Art No. 5914520




Art No. 5914090




Art No. 5914510




Art No. 5914912




Art No. 5914200



Art No. 5914170




Art No. 5914185




Art No. 5914190

NEW




Art No. 5914195

NEW




Art No. 5914122

NEW



Art No. 5914522

NEW



Art No. 5914822



Art No. 5913650

1/2.7"
1920 × 1080 pixels

1:1

Art No. 5914321

1/2.8"
2048 × 1536 pixels

1:1

Art No. 5914520

1/2.5"
2529 × 1944 pixels

1:1

Art No. 5914090

1"
4096 × 2160 pixels

1:1

Art No. 5914510

2/3"
2448 × 2048 pixels

1:1

Art No. 5914912

1/1.7"
4000 × 3000 pixels

1:1

Art No. 5914200

1"
5440 × 3648 pixels

1:1

Art No. 5914170

1/2.7"
1920 × 1080 pixels

1:1

Art No. 5914185

1/1.8"
1920 × 1080 pixels - Live
2592 × 1944 pixels - Foto

1:1

Art No. 5914190

1/1.9"
1920 × 1080 pixels

1:1

Art No. 5914195

1/1."
3840 × 2160 pixels

1:1

NEW

Art No. 5914122

1/2.8"
1920 × 1080 pixels

1:1

NEW

Art No. 5914522

1/2.8"
2592 × 1944 pixels

1:1

NEW

Art No. 5914822

1/0.25"
2592 × 1944 pixels

1:1

NEW

Art No. 5914250

1/2.3"
4928 × 4928 pixels

1:1

Overview: MikroCams



Technical Data

Art. No.	Product name	Resolution pixels	pixels size µm	Zoll "	Length mm	Width mm	Area mm²	Diagonal mm	Aspect ratio	Mode
5913650	MikroOkular Full HD Okularkamera	1920 x 1080	3.0	1/2.7	5.76	3.24	18.66	6.61	16:9	
5914321	MikroCam SP 3.1 Mikroskopkamera	2048 x 1536	2.5	1/2.8	5.28	3.84	20.28	6.53	4:3	
5914520	MikroCam SP 5.0 Mikroskopkamera	2592 x 1944	2.2	1/2.5	5.70	4.30	24.51	7.14	4:3	
5914090	MikroCamII 9MP 4K 1" Mikroskopkamera	4096 x 2160	3.45	1.0	14.13	7.45	105.27	15.97	~ 17:9	
5914510	MikroCamII 5MP HIS Mikroskopkamera	2448 x 2048	3.45	2/3.0	8.45	7.07	59.74	11.02	5:4	
5914912	MikroCamII Mikroskopkamera 12MP USB 3.0	4000 x 3000	1.85	1/1.7	7.40	5.55	41.07	9.25	4:3	
5914200	MikroCamII 20 MP 1" Mikroskopkamera	5440 x 3648	2.4	1.0	13.06	8.76	114.41	15.73	3:2	
5914170	MikroCam MINI Full HD HDMI Mikroskopkamera	1920 x 1080	2.9	1/2.8	5.57	3.13	17.43	6.39	16:9	
5914185	MikroCam PRO HDMI Mikroskopkamera 5MP	1920 x 1080 2592 x 1944	2.4	1/1.8	6.22	4.67	29.05	7.78	16:9 4:3	Live Photo
5914190	MikroCam PRO HDMI Autofocus Mikroskopkamera	1920 x 1080	3.75	1/1.9	7.20	4.05	29.16	8.26	16:9	
5914195	MikroCam PRO HDMI 4K	3840 x 2160	2.0	1/1.8	7.68	4.32	33.18	8.81	16:9	
5914122	MikroCam SF Full HD	1920 x 1080	2.9	1/2.8	5.57	3.13	17.43	6.39	16:9	
5914522	MikroCam SF 5.0 MP	2592 x 1944	2.0	1/2.8	5.18	3.89	20.15	6.48	4:3	
5914822	MikroCam SF UHD	3840 x 2160	1.62	1/2.5	6.22	3.50	21.77	7.14	16:9	

NEW